System MG

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Co:	~62,0
Cr:	30,0
Mo:	5,5
Si:	1,0
Mn:	0,6
C:	0,6
Sonstige:	< 0,1

CE

0434

Eigenschaften

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System MG is a partial alloy on cobalt base

DIN EN ISO 22674:2007, Type 5

Modelling and Sprue System:

The modellation of the casting objects happens as usual. 2 casting channels with a thickness of 3 - 3,5 mm are normally sufficient for the sprue system. These will be added fluently at the massiest parts of the modellation.

Investing and Preheating:

Any commercial phosphate bonded investments, which are suitable for a preheating temperature of 950°C, are qualified for the investing (please follow the handling instructions of the investment-producer).

Preheating of the muffle according to the instructions of the investment producer to an end temperature of 950°C – 1050°C.

Melting and Casting:

! Use only clean ceramic crucibles.

! Use an extra crucible for every alloy.

! Only use new metal for optimum casting results.

! Do not use melting powder.

Pre-melting of the alloy in the ceramic crucible. The muffle will be placed into the centrifuge, before the last cylinder has collapsed. Cast alloy until the shadow in the center disappears. The alloy is now completely melted and the casting procedure should be started immediately.

! Do not overheat the alloy.

Please follow the instructions of the manufacturer while melting with acetylene/oxygen. A cleanly adjusted flame prevents contamination of the alloy. The casting procedure will be started, as soon as the alloy is melted and moves slightly under flame pressure. After the casting cool the muffle to room temperature and deflasking.

! No water quench.

Soldering:

Any commercial high-melting cobalt-chrome solders may be used for the soldering of the alloys.

If the alloy is handled according to the corresponding instructions, it can be polished and grinded easily. A good high gloss polishing guarantees an optimum density and best biocompatibility.

Guarantee: All recommendations regarding the handling are based on our own experiences and Therefore they may only be taken as guide values. Our products are subject to a continuous advancement. Construction and composition are subject to alteration.