# System Duro

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
## Non-precious dental alloy on cobalt/chrome base

### Component of the alloy in %

Co:	Rest
Cr:	25,0
Mo:	1,0
Si:	1,0
C:	0,04
W:	9,0
Fe:	1,0
Nb:	2,0
V:	2,0
Others:	< 0,1

### Characteristics

Hardness after casting	> 340 HV 10
Specific gravity	8,4 g/cm <sup>3</sup>
Solidus/Liquidus	1.370-1.385°C
Casting temperature	1.445°C
WAK (20°C-600°C)	14,1 µm/mK
Delivery form	Cylinder ø 8 x 15 mm
Packaging unit	2,2 lb (1000 g)

System Duro is a nickel- and beryllium-free cobalt/chrome alloy

DIN EN ISO 22674:2006, Type 5

### Modelling and Sprue System

To guarantee a safe flow of the modelling, the thickness of the cap should not fall below 0,3 - 0,4 mm. The sprue system of the objects happens as usual with sprue reservoir or crossbar (crossbar should be divided for large bridges). Massive bridge parts should be provided with air outlets.

### Investing and Preheating

Any commercial phosphate bonded investments, which are suitable for a preheating temperature of 950 °C, are qualified for the investing (please follow the handling instructions of the investment-producer).

Recommendation: Large casting objects should only be heated linear!

### Casting

The optimum preheating temperature of the muffle is 900°C. A ceramic crucible must be used for the melting of the alloy (no graphite crucible!). Casting steam should be exhausted.

### Melting of the alloy

**with open flame (acetylene/oxygen or propane/oxygen)** with neutral flame without melting powder. When the casting cubes are completely melted, (melting moves under flame pressure), the casting procedure should be started immediately. The oxide layer may not crack, because otherwise alloy components might burn, and this leads to changes of the characteristics of the alloy!

When melting the alloy **in the high frequency centrifuge**, the casting procedure will be started immediately after all cubes have collapsed and as soon as the shadow over the melting is disappeared. Even here the oxide layer may not crack!

### Cool muffle slowly in the air.

The elaboration of the modelling happens with suitable milling cutters and Al-Oxyd stones. The thickness of the cap should not go below 0,2 - 0,3mm after finalization.

### Recommended for single use

You should always use new casting material, as after repeated casting of the alloy important elements will evaporate and an optimum bonding between metal and ceramic is no longer guaranteed.

### Firing of the ceramic

An oxide firing is not necessary. If an oxide firing is desired, then under atmosphere -10 min at 980°C (wash brand is not applicable for bonder brand). Then sandblast with aluminum oxide 110µm. Never put non-precious dental alloys in acid-stain bath. Long-term-cooling used for all firings except opaque.

### Soldering and Welding

Commercially available solder on cobalt base is suitable for soldering. Never use gold- or palladium solder for the soldering of NE objects among themselves

**If the alloy is handled according to the corresponding instructions, it can be polished and grinded easily. A good high gloss polishing guarantees an optimum density and best biocompatibility.**

### Cleaning

in ultrasonic bath or with steam blasting.

### Guarantee

All recommendations regarding the handling are based on our own experiences and tests. Therefore they may only be taken as guide values. Our products are subject to a continuous advancement. Construction and composition are subject to alteration.

### Security advice

Metal dusts are harmful to your health. Therefore suitable aspiration and a dust mask are necessary when finalizing and blast cleaning.

Our information and recommendations are subject to the current state of the scientific and technical knowledge and are to our knowledge and experiences presently to be seen as correct. They do not include any attribute confirmations. The forthcoming version replaces all former specifications.