System KN

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OE

0434

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Ni :	-62,0
Cr:	25,0
Mo:	11,0
Si:	1,5
Mn:	<0,1

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Non-precious dental alloy on nickel base

OE

0434

Component of the alloy in %:

Ni:	-62,0
Cr:	25,0
Mo:	11,0
Si:	1,5
Mn:	<0,1

Hardness after casting	190 HV 10
Hardness after firing	185 HV 10
Specific gravity	8,2 g/cm ³
Solidus/Liquidus	1.260-1.350°C
Casting temperature	1.420°C
WAK (20°C-600°C)	14,0 µm/mK
Delivery form	cylinder ø 9,5 x 11 mm
Packaging unit	2,2 lb (1000 g)

System KN is a beryllium free dental alloy on chrome/nickel base

Modelling and Sprue System:

To guarantee a safe flowing of the modelling, the cap thickness should not fall below 0,3 - 0,4 mm.

The sprue system of the objects happens as usual with sprue reservoirs or crossbar (crossbar should be divided for large bridges). Massive bridge parts should be provided with air outlets.

Investing and Preheating:

Any commercial phosphate bonded investments, which are suitable for a preheating temperature of 1050°C, are qualified for the investing (please follow the handling instructions of the investment-producer).

Preheating of the muffle according to the instructions of the investment producer to an end temperature of 900°C – 1050°. Large casting objects should only be heated linear (no speed-investment)!

Melting and Casting:

! Use only clean ceramic crucibles.

! Use an extra crucible for every alloy.

! Only use new metal for optimum casting results. After repeated casting of the alloy important elements could evaporate and an optimum bonding between metal and ceramic is not longer guaranteed.

! Do not use melting powder.

Pre-melting of the alloy in ceramic crucible.

The muffle will be placed into the centrifuge, before the last cylinder has collapsed.

Cast alloy until the shadow in the center disappears.

The alloy is now completely melted and the casting procedure should be started immediately.

! The oxide skin may not crack (overheating)

Please follow the instructions of the manufacturer while melting with acetylene/oxygen.

A cleanly adjusted flame prevents contamination of the alloy

The casting procedure will be started, as soon as the alloy is melted and moves slightly under flame pressure.

After the casting cool the muffle to room temperature and deflasking.

! No water quench.

The grinding of the modelling happens with suitable milling cutters, stones and sinter diamonds. The thickness of the cap should be at least 0,2 - 0,3mm

Firing of the ceramic

An oxide firing is not necessary. If an oxide firing is desired, then under atmosphere (wash brand is not applicable for bonder brand). Then sandblast with aluminum oxide 110µm and clean with distilled water.

Please note WAK-values NE/ceramic!

Longtime cooling is not necessary, however it is recommended to remove the work from the tray after the firing procedure after the red heat is disappeared.

Soldering and Welding

Commercially available solder on cobalt base is suitable for soldering.

Never use gold- or palladium solder for the soldering of NE objects among themselves

Cleaning in ultrasonic bath or with steam blasting.

Guarantee: All recommendations regarding the handling are based on our own experiences and tests. Therefore they may only be taken as guide values. Our products are subject to a continuous advancement. Construction and composition are subject to alteration.