System NE

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Non-precious dental alloy on cobalt/chrome base

CE

0434

Component of the alloy

Co:	~65,0
Cr:	20,0
W:	6,5
Mo:	6,5
Si:	0,85
Fe:	< 0,5
Mn:	0,85

Characteristics

Specific gravity	8,4 g/cm ³
Solidus/Liquidus	1.280-1.350 °C
Casting temperature	1.420 °C
WAK (20°C-600°C)	14,8 µm/mK
Delivery form	cylinder ø 9,5 x 11 mm
Packaging unit	2,2 lb (1000 g)

System NE is a nickel- and beryllium-free cobalt/chrome alloy without carbon

Modelling and Sprue System:

To guarantee a safe flow of the modelling, the thickness of the cap should not fall below 0,3 - 0,4 mm.

The sprue system of the objects happens as usual with sprue reservoirs or crossbar (crossbar should be divided for large bridges). Massive bridge parts should be provided with air outlets.

Investing and Preheating:

Any commercial phosphate bonded investments, which are suitable for a preheating temperature of 1050°C, are qualified for the investing (please follow the handling instructions of the investment-producer).

Large casting objects should only be heated linear (no speed investment)!

Casting

The optimum preheating temperature of the muffel is 950°C.

A **ceramic crucible** must be used for the melting of the alloy (**no graphite crucible!**). Casting steam should be exhausted.

Melting of the alloy with open flame (acetylene/oxygen or propane/oxygen) with neutral flame without melting powder. When the casting cubes are completely melted, (melting moves under flame pressure), the casting procedure should be started immediately. The oxide layer may not crack, because otherwise alloy components might burn, and die change the characteristics of the alloy!

When melting the alloy in the high frequency centrifuge, the casting procedure will be started immediately after all cubes have collapsed and as soon as the shadow over the melting is disappeared. Even here the oxide layer may not crack!

Cool muffle slowly in the air.

You should always use new casting cubes, as after repeated casting of the alloy important adhesive oxide former will evaporate and an optimum bonding between metal and ceramic is not longer guaranteed.

The elaboration of the modelling happens with suitable milling cutters, stones and diamonds. The thickness of the cap should be at least 0,2 - 0,3mm.

Firing of the ceramic

An oxide firing is not necessary. If an oxide firing is desired, then under atmosphere (wash brand is not applicable for bonder brand). Then sandblast with aluminum oxide 110µm and clean with distilled water.

Soldering and Welding

Commercially available solder on cobalt base is suitable for soldering.

Never use gold- or palladium solder for the soldering of NE objects among themselves

Cleaning in ultrasonic bath or with steam blasting.

Guarantee: All recommendations regarding the handling are based on our own experiences and tests. Therefore they may only be taken as guide values. Our products are subject to a continuous advancement. Construction and composition are subject to alteration.